Adult Resuscitation Trolley Checklist, with additional paediatric BLS equipment.

	Mon	Tue	Wed	Thurs	Fri	Sat	Sun
Date:							
Time:							
Print name:							
Signature:							
Top/Side of trolley							
Top Defibrillator with 1xAdult/Child pads- passed daily self-check	Run Op Check						
Top ECG dot electrodes (1x pack of 3)							
<i>Top</i> Sharps Bin							
Left side CD sized oxygen cylinder (full)							
Left side Bougie Unopened							
Left side Suction Catheters size 12Ch & 14Ch							
Left side Gloves (Medium and Large/XL)							
Right side Suction Unit (Plugged in and tested)							
Drugs Drawer		The Princess Alexandra Hospital NHS					
In-date Cardiac Arrest drugs as per trolley list			76.11			NHS T	
Disposable razor			200: EEL -				
Airway Drawer		37 %	39 0 966	F D	:4 . 4: 7	Facilian Daa	
Sealed (Blue) Airway Tray & in-date	For Resuscitation Trolley Resources						
Breathing Drawer			.5175	< Scan	the QR Co	nda	
Disposable Adult Bag-Valve-Mask & Disposable Paediatric Bag- Valve-Mask (500ml)				← Scarr		Jue	
Adult Mapelson-C Breathing Circuit							
Adult & Paediatric Non-Rebreathe Mask]	14/	1/1 1:			
Circulation/IV Access drawer			ward	d/Location:		<u></u>	
Sealed (Red) IV access tray and in-date			Defibrillato	r: Heartstart Intr	epid/Philips X	L+/CR2 (Circle)	
IV Fluids Drawer		Heartstart Intrepid Defibrillator: Additional Operation check required weekly; (guidance available via QR Code above) Please perform operation check every Monday					
2x blood & 2x IV administration giving sets							
2x 1000ml bags of 0.9% saline							
2x 500ml bags of colloid solution							
1 x Pressure bag			<u> </u>	ease perioriii operati	on check every w	<u>onday</u>	
Miscellaneous Drawer				Heartstart Intrep		ıy	
Suction catheters - 2 x size 12 & 14]		Operation Chec	k		
Spare suction unit liner, tubing & yankeur]		Tiek When Com	olotod:		
1x spare set of compatible Adult/Child defibrillator pads]		Tick When Comp	pieteu.		
Spare ECG dot electrodes (4x pack of 3)]	L				
Ligature Cutter]					
Spare trolley seals]					
Trolley <u>cleaned, dusted</u> & <u>sealed</u> after check completed	Seal number:	Seal number:	Seal number:	Seal number:	Seal number:	Seal number:	Seal number:
If re-check required please add seal number and initial							

Adult Resuscitation Trolley Checklist, with additional paediatric BLS equipment.