

Seal Your Trolley...

Your Avalo Resus Trolley uses our **Avalo Breakaway Lock Seal** in Red - Item No. **7923** (available in **6 different colours**). Below is a guide recommending correct practice in order to **keep your Trolley tamper evident**.



1. Close all drawers and lock your cart by pushing the blue handle inwards and then upwards. You will see a small hole for which to insert your seal.



2. Thread your seal through the hole visible under the cart handle.



3. Loop the tail back through the security seal and record its number. Recording a Seal's number and logging in a check book saves valuable time checking each draw every day.


Avalo Breakaway Lock Seal

Pull-Fly Baby Seals fit small boxes, Resus/Crash Trolleys and any area with limited space, where other seals may have difficulty sealing or breaking.

Our Baby Seals are the only seal approved for use with the Avalo breakaway locking handle.

Breaking Strength: **LOW**
100 Seals Per Pack

£19.27
per pack



How to order...

1. Make sure you have an official **Purchase Order** number from your supplies department
2. Make sure you have the correct **Item Number** corresponding with the colour you require
3. **Call** 0800 525 076 or **email** your order to info@distinctivemedicalproducts.com or **Fax** your order to 0800 590 014



Call Today: 0800 525 076
NEW Website & Brochure available Jan 2013

distinctive
medical